

PURCHASE ORDER
PROVINCIAL GOVERNMENT OF CAMARINES NORTE

Supplier :	<u>SOPHIES GENERAL MERCHANDISE</u>	P.O. No. :	<u>25102689</u>
Address :	<u>P-5, Bautista, Labo, Camarines Norte</u>	Date :	<u>12-11-25</u>
Telephone No. :	<u>0968-413-2568</u>	Mode of Procurement :	<u>SVP</u>
TIN :	<u>184-646-631-000</u>		

Gentlemen: Please furnish this office the following articles subject to the terms and conditions contained herein:

Place of Delivery: <u>LDH</u>	Delivery Term: <u>FOP</u>
Date of Delivery: <u>12/19/25</u>	Payment Term: <u>FULL</u>

Item No.	Quantity	Unit Issue	ITEM DESCRIPTION	Amount	
				Unit Cost	Total Cost
1	1	unit	PATIENT MONITOR Technical Specification ECG: 3/5 Leadmodes, 2-6 channel display HR: 15-300 bpm, accuracy ± 1 bpm or ±% Gain: 2.5, 5, 10, 20 mm/mV RESP (Respiration) Method: Impedance (RA-LL) RateL: 0-150 rpm, accuracy ±2 rpm SPO2: Range: 0-100% rpm, accuracy ±2% (70-100%) Pulse Rate: 20-250 bpm NIBP (Non-invasive BP) Mode: Manual, Auto, STAT Range: Adult 10-270 mmHg, Pediatric 10-200 mmHg, neonatal 10-135 mmHg Accuracy: ±5 mmHg TEMPL Channel: 2, range 0-50°C, accuracy ±0.1 °C RECORDER: 2 Channel, real-time/event/manual modes Paper Speed: 25-50mm/s 12-CHANNEL ECG (Cardio Fax M 1350A) Key Features: Quality Diagnosis- Advance filtering reduces noise; ECAPS 12C program provides 200+ findings for accurate faster diagnosis. Quality Examination- 8" FTF display for 12-lead monitoring; accurate measurements; stores up to 3 minutes CG waveforms Comfortable Operation- Simple design with user-friendly interface. Identify Invisible Ischemia- Synthesized ECG adds 6 leads + ST-chart for right ventricular & posterior wall analysis without extra procedures Digital Workspace- polaris.one platform enables storage sharing and access of patient data across departments with HIS integration. Quality Accessibility-Data output via LAN (wired/wireless), DICOM compatibility with RIS/PACS <i>for Capital Outlay/Medical Equipment use</i>	64,000.00	64,000.00
2	1	unit	12-CHANNEL ECG (Cardio Fax M 1350A) Key Features: Quality Diagnosis- Advance filtering reduces noise; ECAPS 12C program provides 200+ findings for accurate faster diagnosis. Quality Examination- 8" FTF display for 12-lead monitoring; accurate measurements; stores up to 3 minutes CG waveforms Comfortable Operation- Simple design with user-friendly interface. Identify Invisible Ischemia- Synthesized ECG adds 6 leads + ST-chart for right ventricular & posterior wall analysis without extra procedures Digital Workspace- polaris.one platform enables storage sharing and access of patient data across departments with HIS integration. Quality Accessibility-Data output via LAN (wired/wireless), DICOM compatibility with RIS/PACS <i>for Capital Outlay/Medical Equipment use</i>	355,000.00	355,000.00


Total amount in words: **Four Hundred Nineteen Thousand Pesos Only** 419,000.00

In case of failure to deliver within the time specified above, penalty of one-tenth (1/10) of one (1) percent of every day of delay shall be imposed.

Very truly yours,

ENGR. JOSEPH V. ASCUTIA
Acting Governor

Conforme:


SOFIA P. CALMATEO
 Signature over printed name of Supplier
12-12-25
 Date